

	SCALE IN	DIMENSIONS A		
	INCHES OF WATER OR			MERCURY REQ'D
MODEL	MERCURY	W/M	D	TO FILL (WT.)
1235-8	0-8	17 3/8	18 1/2	12 OZ
1235-20	0-20	29 5/16	33 9/16	18 OZ
1235-36	0-36	45 1/16	51 15/16	26 OZ

©D = CRITICAL DIMENSION
STANDARD TOLERANCES UNLESS NOTED:
ALL DECIMAL DIMENSIONS ± .005
ALL ANGLES ± 1'

			DATE	NAME		MATERIAL			
			7-24-01						
			DWN BY	1235 MANOMETER		FINISH			
			CHKD						
						DWYER INSTRUMENTS, INC.			
NO.	CHANGES	BY/DATE	APPD		10100000	MICHIGAN CITY, INDIANA 46360 U.S.A.			
		ACAD2002	FR. NO.						
	NOTICE: This drawing and the principles and elements of design embodied therein are the exclusive property of DWYER INSTRUMENTS, INC. and are not to be communicated, disclosed, reproduced or used except as previously authorized in writing by such corporation and must not be submitted to outside parties for examination writing thou the writing corporation.								